



# Designing for Commercial Scale: A Rapid and Scalable Lentiviral Manufacturing Approach

Cédric Rousseaux,

Director Innovation, Analytical and Process Development

March 2026

The information contained herein is considered to be confidential information and should be handled accordingly. For planning purposes only. © 2026 Proprietary and Confidential Information of SK pharmteco Inc. and affiliates. | All rights reserved.

# Built on Stability, Committed to Your Success

Investment to support our clients' needs

CDMO built to deliver from  
early development  
through commercialization  
across a global network

**Viral Vectors** | Small Molecules  
*AAV, LVV, Adeno*



# Strategic Facilities Mitigating Potential Supply Chain Risks





Committed to support our clients **deliver life-impacting medicines for patients** around the world

## Mission

We are committed to changing the world through our unwavering dedication to lifesaving therapies. We enable our partners with cutting-edge technologies, expertise, and customer centricity, allowing them to bring breakthrough treatments to patients and deliver happiness and health to stakeholders everywhere

## Vision

The most trusted global partner in the delivery of innovative medicine

# Manufacturing Capability Snapshot

**3** Plug-and-play viral vector manufacturing platforms

Specializing in **AAV, LV-based, Adeno,** and other viral modalities

LVV *ex vivo* to *in vivo* approach



**2** Locations:  
France and US



**100+**  
Number of total cGMP batches produced over AAV, Lentivirus, and Adenovirus



**9** Vector suites +  
Formulation and  
Fill/Finish suites



**3-5** Week Avg.  
Batch Production



**8** Months from plasmid to cGMP release using our LVV and AAV platforms



Up to: **500 L**  
Adherent process

**2 x 1,000 L**  
Suspension process

# Overview & Technical Narrative

From Proven *Ex vivo*  
Manufacturing to *In vivo*  
Drug-Product Readiness

1

## Defining Our LentiSure™ Platform:

A proven and scalable *ex vivo* LVV manufacturing platform built on years of clinical experience.

2

## Our Approach from *Ex vivo* to *In vivo*:

Building on our *ex vivo* LVV track record, we are evolving to drug-product standards with strengthened purity, potency, and CMC control.

3

## Redefining LVV Downstream for *In vivo* Applications:

Addressing yield limitations and impurity constraints to meet drug-product-level expectations.

4

## Scalable Performance & Manufacturing Readiness:

Demonstrating robustness, cost efficiency, and yield from bench scale through 50 L.

From Proven  
*Ex vivo*  
Manufacturing  
to *In vivo*  
**Drug-Product  
Readiness**

1

**Defining Our LentiSure™ Platform:**

A proven and scalable *ex vivo* LVV manufacturing platform built on years of clinical experience.

2

**Our Approach from *Ex vivo* to *In vivo*:**

Building on our *ex vivo* LVV track record, we are evolving to drug-product standards with strengthened purity, potency, and CMC control.

3

**Redefining LVV Downstream for *In vivo* Applications:**

Addressing yield limitations and impurity constraints to meet drug-product-level expectations.

4

**Scalable Performance & Manufacturing Readiness:**

Demonstrating robustness, cost efficiency, and yield from bench scale through 50 L.

# USP: Engineering Biological Control for Commercial Scalability

A modular, fully single-use upstream platform designed for robustness, flexibility, and scale consistency



## Controlled Biological Foundation

- In-house cGMP HEK293T MCB / WCB
- Serum-free culture
- Standardized quadruple-plasmid system
- Proprietary transfection reagent



## Robust & Modular by Design

- DoE-optimized DNA/cell and plasmid ratios
- Standardized transfection parameters
- Fully single-use closed architecture
- Compatibility with customer cell lines and plasmids




## Scalable Performance Demonstrated

- Consistent titers across scales
- Harvest titers up to  $5.9 \times 10^7$  IG/mL
- No loss of infectivity during scale-up
- Demonstrated scalability to 200 L


# Productivity Across Volumes and Transgenes

LentiSure maintains infectious titer and productivity from bench to pilot scale

GFP lentivirus vector scale up



Lentivirus vector productivity at 200 L scale



## Demonstrated scalability:

- Comparable infectious titers from 250 mL to 200 L
- No significant loss of productivity during scale-up
- Stable p24 levels across volumes

## Robustness Across Transgenes:


- Consistent performance across multiple constructs
- Titers up to  $5.9 \times 10^7$  IG/mL at 200 L
- No construct-specific productivity collapse

# DSP: High-Recovery, Scalable Purification by Design

Single-use DSP strategy engineered to preserve infectivity and maximize manufacturing efficiency

## 1 Protecting Vector Integrity Early

- Optimized single-use clarification train
- High-capacity filtration strategy
- 0.2  $\mu\text{m}$  filtration to protect AEX membrane
- Controlled shear exposure




## 2 Membrane-Based AEX for Scalable Recovery

- Single-use membrane chromatography
- High load capacity
- Reduced residence time, minimizing vector degradation
- High recovery of infectious particles



## 3 Efficient Concentration & Manufacturing Readiness

- Single-use TFF/Diafiltration
- Controlled buffer exchange
- Concentration without infectivity loss
- Removal of residual DNA and process-related impurities
- Reduced process time and footprint



# Maintained Infectivity and Controlled Recovery Through DSP

From harvest to drug product without significant infectivity loss



Scan to explore a real-world LentiSure case study.

Follow us on LinkedIn for our upcoming white paper on ex vivo LVV manufacturing.

## Harvest Performance

Comparable harvest titers:


- $3.7 \times 10^7$  IG/mL at 50 L
- $2.9 \times 10^7$  IG/mL at 200 L
- No major scale-related drop

## Drug Substance Intermediate Concentration

Effective downstream concentration

- Infectious titer maintained post-purification:
  - $4.3 \times 10^9$  IG/mL at 50 L
  - $4.0 \times 10^9$  IG/mL at 200 L

## Process Efficiency



- Controlled recovery from harvest to DS intermediate
- Stable across scale

# From Proven *Ex vivo* Manufacturing to *In vivo* **Drug-Product Readiness**

1

## Defining Our LentiSure™ Platform:

A proven and scalable *ex vivo* LVV manufacturing platform built on years of clinical experience.

2

## Our Approach from *Ex vivo* to *In vivo*:

Building on our *ex vivo* LVV track record, we are evolving to drug-product standards with strengthened purity, potency, and CMC control.

3

## Redefining LVV Downstream for *In vivo* Applications:

Addressing yield limitations and impurity constraints to meet drug-product-level expectations.


4

## Scalable Performance & Manufacturing Readiness:

Demonstrating robustness, cost efficiency, and yield from bench scale through 50 L.

# Shift in Regulatory Burden From Drug Substance to Drug Product

Implications for Manufacturing, CMC, & Analytical Strategy



## Key Milestone:

SK pharmteco successfully manufactured and released both **ex vivo** and **in vivo** lentiviral products for clinical use, demonstrating capability across drug substance and drug product manufacturing.

# In vivo CAR-Ts: Barriers to Entry

## Identifying the Key Biological and Manufacturing Bottlenecks

### Tropism & Cell Selectivity

- Off-target transduction of non-T-cell populations (e.g., macrophages or NK cells) may trigger uncontrolled cytokine release and systemic inflammation.
- Unintended antigen recognition on healthy tissues (e.g., low-level HER2 expression) can result in fatal on-target/off-tumor toxicity.
- Achieving strict T-cell-specific tropism is essential to minimize immune-related adverse events.

### Regulating Expression

- Current LVVs often rely on strong, constitutive promoters, leading to poorly regulated expression kinetics.
- I. Underexpression → loss of potency & persistence
  - II. Overexpression → cytokine release syndrome or tissue damage
  - III. Variable expression → heterogeneous potency & unpredictable safety profiles.

### Manufacturing & Dosing Challenges

- *In vivo* LVVs represent the final drug product (DP) rather than an intermediate, imposing stricter CMC and analytical expectations.
- Requires robust vector quality and durability to ensure consistent *in vivo* performance.
- Variability in transduction efficiency, T-cell expansion kinetics, and host immunity affects potency and comparability across patients. Neutralization of transduced cells before expansion.

# From Proven *Ex vivo* Manufacturing to *In vivo* **Drug-Product Readiness**

1

## Defining Our LentiSure™ Platform:

A proven and scalable *ex vivo* LVV manufacturing platform built on years of clinical experience.

2

## Our Approach from *Ex vivo* to *In vivo*:

Building on our *ex vivo* LVV track record, we are evolving to drug-product standards with strengthened purity, potency, and CMC control.

3

## Redefining LVV Downstream for *In vivo* Applications:

Addressing yield limitations and impurity constraints to meet drug-product-level expectations.

4

## Scalable Performance & Manufacturing Readiness:


Demonstrating robustness, cost efficiency, and yield from bench scale through 50 L.

# Dissecting AEX Chromatography

## Why Traditional AEX Constrains LVV Manufacturing

**Goal:** Identify why resin-based AEX limits LVV recovery and scalability and what must change to unlock higher yields.

### Traditional IEX Process



### Comparison of IEX Technologies

Attribute	Packed-Bed Resin	Membrane Adsorber
Mass Transfer	Diffusion-limited	Convection-driven
Recovery	15–40%	60–90%
Throughput	Slow, batch-limited	>5× faster
Back Pressure	Increases at High load	Minimal
Residence Time	2-4 min	0.1 min
Cost & Flexibility	High COGs, rigid hardware	Disposable, cost-efficient

### Impact of Process Optimization on LVV Elution

Elution	VP/mL	Total VP	%VP of Load
Traditional	1.1E+08	4.6E+08	17%
Optimized	4.6E+08	1.9E+09	70%


**Even with optimization, resin-based AEX remains diffusion-limited thus capping LVV recovery, extending cycle times, and driving high COGs.**

# Tripling LVV Recovery with Optimized IEX Chromatography


Driving Recovery from 20% to 70% Through High-Throughput Convective IEX Membrane Chromatography

**Goal:** To overcome diffusion-limited recovery, we evaluated convective membrane-based IEX as an alternative capture strategy.


## SK pharmteco's Improved Process:




## Improved Yields



## Reduced Residuals



## Intact Virus



## Summary:

- Achieves 98% HCP clearance (99.7% from harvest).
- **> 85% TU recovery** across runs.
- **≥ 70% overall LVV yield.**
- **<2h processing time**, enabling simplified operations.

**Convective membrane flow eliminates diffusion constraints, enabling high-throughput viral capture with improved impurity resolution.**

# From Proven *Ex vivo* Manufacturing to *In vivo* **Drug-Product Readiness**

1

## **Defining Our LentiSure™ Platform:**

A proven and scalable *ex vivo* LVV manufacturing platform built on years of clinical experience.

2

## **Our Approach from *Ex vivo* to *In vivo*:**

Building on our *ex vivo* LVV track record, we are evolving to drug-product standards with strengthened purity, potency, and CMC control.

3

## **Redefining LVV Downstream for *In vivo* Applications:**

Addressing yield limitations and impurity constraints to meet drug-product-level expectations.






4

## **Scalable Performance & Manufacturing Readiness:**

Demonstrating robustness, cost efficiency, and yield from bench scale through 50 L.

# Consistent LVV Performance from Bench to 50 L Scale

Optimized IEX Processing Delivers Robust Recovery, Purity, and Fast Processing Across Scales



Data shown represent representative internal process development runs; performance may vary by vector construct and upstream conditions.

**Total DS Processing Time from Clarification to SF: 6 – 8 hrs**

**Optimized IEX processing maintains consistent LVV recovery, impurity clearance, and processing times from benchtop to 50 L, demonstrating scalability and manufacturing robustness.**

# Go Deeper: *In Vivo* & *Ex Vivo* LVV Manufacturing

Explore the full technical strategy behind our optimized LVV platform

Prefer to read?

**Download our white paper:**

*Revolutionizing Lentiviral Vector Manufacturing for In Vivo & Ex Vivo CAR-T*

Scan to access



Prefer to listen?


**Watch our on-demand webinar:**

*Enhancing LVV Manufacturing: Eliminating the Barriers to In Vivo and Ex Vivo CAR-T Delivery*

Scan to watch



# From Platform Excellence to *In Vivo* Drug-Product Readiness





# Thank You

**Meet us at booth #203** to learn how we can support your program requirements.



[www.skpharmteco.com](http://www.skpharmteco.com)




LinkedIn [@SK pharmteco](#)



# Lentiviral Vector (LVV) Downstream Is Yield-Limited by IEX Recovery

IEX remains the dominant loss point despite downstream optimization



**Even with downstream optimization, IEX recovery remains the dominant constraint on LVV yield (typically 15-30% LVV recovery).**